


Preparation of all-trans retinal food

Jelly hubertus Maria Soffers

 0000-0002-1051-7375

j.soffers@umkc.edu

1 Retinal Food Preparation

Protocol modified from:

- Optogenetics teaching protocol
- Molasses food recipe

2 Preparation of All-Trans Retinal (ATR)

- Suspend ATR in ethanol to make a 250 mM stock solution.
- Prepare fly food with a 1:200 dilution → final concentration 0.2 mM.
- ATR does not fully dissolve, resuspend as much as possible.
- Store 1 mL aliquots at -20°C (stable indefinitely).

Light Sensitivity of ATR

All-trans retinal is highly light sensitive. Avoid exposure to strong white light and keep solutions covered with foil whenever possible.

3 Fly Media Options

3.1 Option 1: Pre-made fly media

- Purchase fly media online Archon Scientific
- Supplement with ATR as needed
- **Clark lab:**
 - Drip ATR onto food surface
- **Lacin lab:**
 - Scrape food from ~40 vials (~5 mL each)
 - Microwave to melt and cool
 - Add 2 vials of 250 mM ATR to 40 vials

3.2 Option 2: Prepare your own fly media

Consumables and Equipment:

- Please refer to the vendors listed by the Bloomington Drosophila Stock Center for supplies and equipment. Vendors Fly supplies.
 - Genesee Scientific is amongst the most competitively priced vendors for fly vials, plugs and food.
- A simple kettle on heat plate or turkey roaster oven (~\$75.00) are economic solutions for cooking the fly food.
 - Turkey Roaster Oven

Recommended mix:

- Nutri-Fly™ Molasses formulation Genesee Scientific] () (Costs < \$3.00 per liter)
 - Nutru-Fly™
- Or use the recipe below:

3.2.0.1 Yield

- Makes ~1.2 L
- Enough for ~100 vials

3.2.0.2 Ingredients

Ingredient	Amount
Water	1041.5 mL
Yeast Red Star Active Dry)	26.4 g
Yellow cornmeal	40.8 g
Agar	9.6 g
Molasses	64.83 mL
Tegosept (20% in ethanol)	12 mL
Propionic acid	6.0 mL
Extra water	141 mL

Ingerient vendor: We recommend Genesee Scientific
Instructions

1. Add $\frac{1}{4}$ of the water and stir with:
 - agar
 - yeast
 - cornmeal

- molasses (≥ 10 minutes)
2. Add remaining water and boil for 20 minutes, stirring to dissolve lumps
 3. Add remaining cold water to cool faster
 4. When mixture is $< 80^{\circ}\text{C}$, add:
 - acid
 - tegosept
 5. Add 6 mL of 250 mM ATR

☒ Protect from Light During Preparation

After adding ATR, minimize light exposure. Work quickly and keep containers covered with foil when possible.

6. Dispense ~ 5 mL per vial
7. Dry overnight:
 - cover with cheesecloth
 - seal with Press'n Seal
 - cover tray with foil (light protection)
8. Storage: Best used within 4–6 weeks

☒ Storage Tip

Keep prepared vials covered with aluminum foil and store in low-light conditions to preserve ATR activity.